Advanced
Solutions
in Acute
Respiratory
Care

CPAPValve


A complete line of fixed Isobaric CPAP Valves.


Scaradyne

A broad selection of Isobaric CPAP valves allows the clinician to titrate pressure accurately in accordance with the patient's condition.

When used with an integrated, portable CPAP system such as the WhisperFlow generator, you can reduce the patient's work of breathing; improve their O₂ saturation; be pro-active in treating pulmonary edema; prevent the onset of post-operative atelectasis; aid in treatment of ARDS; assists in alleviating tachypnea and aid in weaning off a ventilator.

- Fixed pressure at any flow
- Available in 22mm and 30mm connections
- Sizes: 2.5, 5.0, 7.5, 10.0, 12.5, 15.0, and 20.0cm of H₂O
- Sizes are color coded for easy identification
- Interfaces directly with portable CPAP systems, manual resuscitators and ventilator circuits
- · No adjusting necessary





CPAP Valve

Fixed pressure Isobaric CPAP valves are accurate devices that deliver a constant pressure under any flow rate through the use of precision spring loaded valves.



The family of Isobaric CPAP valves is designed for use with manual resuscitators, portable ventilators and CPAP systems, like the Caradyne WhisperFlow. Pre-assembled and ready to use, these CPAP valves come in seven fixed sizes.

Since adjustable PEEP valves constantly vary as the delivered flow changes, the only way to insure that proper pressure is maintained is to use a fixed "non-flow dependent" CPAP valve. Isobaric CPAP valves are essential when precise pressure is required.

- Clear for easy visualization of moving parts
- · Accurate in all flow conditions
- Full range of sizes for various respiratory conditions
- · Quiet for patient and caregiver
- ISO standard 22mm and 30mm connections
- Easy to use

Product Specifications	
Valve Fittings:	ISO (22mm M) and (30mm F)
Material:	Polypropylene

Product Descriptions				
Part No.:	Description/Flow:	Connection:	Qty/Box:	
8-208503-00	2.5cm H ₂ O	22mm	10	
8-208505-00	5.0cm H ₂ O	22mm	10	
8-208507-00	7.5cm H ₂ O	22mm	10	
8-208510-00	10.0cm H ₂ O	22mm	10	
8-208512-00	12.5cm H ₂ O	22mm	10	
8-208515-00	15.0cm H ₂ O	22mm	10	
8-208520-00	20.0cm H ₂ O	22mm	10	
8-209503-00	2.5cm H ₂ O	30mm	10	
8-209505-00	5.0cm H ₂ O	30mm	10	
8-209507-00	7.5cm H ₂ O	30mm	10	
8-209510-00	10.0cm H ₂ O	30mm	10	
8-209512-00	12.5cm H ₂ O	30mm	10	
8-209515-00	15.0cm H ₂ O	30mm	10	
8-209520-00	20.0cm H ₂ O	30mm	10	
			·	

Caution

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.





Manufactured by:
Caradyne
Parkmore Business Centre
Parkmore West, Galway, Ireland
T: + 353 (0) 91 709010
F: + 353 (0) 91 75829
Email: info@caradyne.com
Web: www.caradyne.com

Caradyne 2055 Executive Drive Indianapolis, IN 46241 USA T: 317.246.9010